




# MEZZO 5 *plus*

## Dynamic Mattress Overlay System

for those patients at **High Risk** or with existing pressure damage.

Together we **will** make an  essential difference

# MEZZO 5 *plus* Mattress

MEZZO 5 *plus* configured with 18 cells which alternate during a 10 minute 2 cell cycle and 4 static head cells. MEZZO 5 *plus* provides effective pressure relief for those patients up to **High Risk** or with existing pressure damage.

## MEZZO 5 *plus* Mattress

- 18 cells operate a clinically proven 10 minute 2 cell alternating cycle.
- Alternating therapy action optimises low interface pressure relief and gentle capillary stimulation.
- 3 Static head cells for patient comfort.



- Timed static function facilitating stability during transfer and patient procedures, auto default to alternating mode after 30 minutes ensures patient safety.
- CPR for rapid mattress deflation during medical emergencies.
- Cable management loops for safety.
- Modular construction reduces maintenance and repair costs.

## MEZZO 5 *plus* Cover

- Multi-stretch to prevent friction & shear.
- Waterproof to prevent fluid ingress.
- Vapour permeable to optimise comfort, preventing skin maceration and breakdown.
- Bacteriostatic, fungistatic and anti-microbial polyurethane fabric, with proven resistance to MRSA.
- High frequency welded seams to prevent fluid ingress.
- Fluid ingress protection flap along the length of the zip.
- Easy to clean - can be wiped down or completely removed and laundered up to 95°C.



# MEZZO 5 *plus* Pump

- Compact & quiet.
- Easy to use.
- Patient weight setting.
- Alternating & Static Function.
- Auto default to alternating mode after 30 minutes ensuring patient safety.
- Function lock facility.
- Audible & visual alarms.
- Alarm mute facility.
- Externally accessed air filter for easier maintenance.



HIGH

High Risk



180kg

Maximum User Weight



2 Year warranty  
Mattress & Cover



High Frequency  
Welded Seams



Engineered For  
Profiling Beds

Specifications			
Mattress:		Pump:	
Dimensions	200 x 88 x 12.5cm (inflated)	Dimensions	29cm x 14cm x 9cm
Weight	3.5 kg	Weight	1.7Kg
Safety Standards	BS 7177: 1996 Crib Source 5	Safety Standards	IEC601/1 BS5724 Part 1
		Voltage:	AC230V 50Hz
		Noise Level:	NC27
		Power Input:	12w
		Fuse:	T1A

## Product Codes:

System EHSM5PS01 • Mattress EHSM5PM01 • Pump EHSM5PP01



**Essential Healthcare Solutions (UK) Ltd**  
Essential Enterprise Village, Field Head Lane,  
Birstall, West Yorkshire WF17 9BN

**Ph** 0845 344 1912  
**Fx** 0845 344 1913

**Em** [info@essential-healthcare.com](mailto:info@essential-healthcare.com)  
**Web** [www.essential-healthcare.com](http://www.essential-healthcare.com)

part of the  **essential** group of companies

